**ThermoFisher**  
SCIENTIFIC

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
Jinsheng Hu

CCS Application Engineer

 The world leader in serving science

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
USP Chromatographic Columns  
<https://www.uspchromcolumns.com/chrom>

SW  
[https://extranet.edqm.eu/publications/Recherches\\_SW.shtml](https://extranet.edqm.eu/publications/Recherches_SW.shtml)

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亲水型反相 Hydrophilic-RP

Accucore C18      Accucore aQ  
 Accucore XL C18      Hypersil GOLD aQ  
                                  Synchronis aQ  
 Hypersil GOLD      Aquasil C18  
 Acclaim 120 C18  
 Synchronis C18      Accucore Polar Premium  
                                  Acclaim PolarAdvantage II  
 Hypersil BDS C18      Acclaim PolarAdvantage  
 Hypersil ODS  
 BetaSil C18      Accucore C30  
 BetaBasic C18      Acclaim C30  
 Umisil C18  
                                  Acclaim Organic Acid  
 C8/C6/C4/C1      Acclaim AmG C18  
                                  Phenyl/PFP

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HILIC

Accucore HILIC  
 Hypersil GOLD Silica  
 Synchronis Silica  
 Hypersil GOLD Amino  
 Synchronis Amino  
 Hypersil APS-2  
 Hypersil GOLD Cyano  
 Hypersil BDS Cyano  
 Hypersil GOLD HILIC  
 Synchronis HILIC  
 Accucore 150-Amide-  
 HILIC  
 Accucore Urea-HILIC  
 Acclaim HILIC-10  
 Betasil Diol

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Hypercarb

Hypercarb

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Mixed-Mode

Acclaim Mixed-Mode  
 WAX-1  
 Acclaim Mixed-Mode  
 WCX-1  
 Acclaim Mixed-Mode  
 HILIC-1  
 Acclaim Trinity P1  
 Acclaim Trinity P2  
 Acclaim Trinity Q1  
 Acclaim Surfactant Plus  
 GlycanPac AXR-1  
**GlycanPac AXH-1**

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Ion exchange

BioBasic SCX  
 Hypersil GOLD SAX  
 Hypersil SAX  
 Hypersil GOLD AX

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HILIC 模式

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Accucore Polar Premium  
Acclaim PolarAdvantage II

Acclaim Organic Acid  
Acclaim AmG C18

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Hypersil GOLD SAX  
Hypersil SAX  
Hypersil GOLD AX  
Acclaim Mixed-Mode WAX-1  
Acclaim Trinity P1/P2

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Accucore Polar Premium  
Acclaim PolarAdvantage II

Hypersil GOLD  
Acclaim AmG C18

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BioBasic SCX  
Acclaim Mixed-Mode WCX-1  
Acclaim Trinity P1/P2

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
## Acclaim Mixed-Mode WAX-1

### Quality Assurance Report

Better Separations Through  
Better Chemistry

**Acclaim™ Mixed-Mode WAX-1**  
5µm 120Å (4.6 X 250 mm)  
Product No. 064985

Date: 01-Oct-21 08:36  
Serial No.: 001973  
Lot No.: 02020135  
Mobile Phase: 50:50 v/v Acetonitrile/50 mM Phosphate buffer, pH 6.0  
Flow Rate: 1.00 mL/min  
Detection: 220 nm  
Injection Volume: 5.0 µL  
Temperature: 30 °C  
Storage Solution: Mobile Phase

No.	Peak Name	Ret.Time (min)	Asymmetry (EP)	Plates/Column (EP)	Concentration (ppm)
1	Uracil	2.78	1.15	23671	100.0
2	Iodide	10.51	1.08	27557	100.0

#### QA Results:

Analyte	Parameter	Specification	Results
Iodide	Efficiency	>=16,811	Passed
Iodide	Asymmetry	0.95-1.32	Passed
Iodide	Pressure	<=1980	938

Production Reference:  
Database: Acclaim7  
Directory: Silica/Silica\_1  
Sequence: 2002544\_064985\_AM  
Sample No: 8  
7.2.10.23925  
065167-02 (QAR)  
Chromatop™ Thermo Fisher Scientific


## Acclaim Mixed-Mode WCX-1

### Quality Assurance Report

Better Separations Through  
Better Chemistry

**Acclaim™ Mixed-Mode WCX-1**  
5µm 120Å (4.6 x 250 mm)  
Product No. 068352

Date: 17-Jun-20 07:15  
Serial No.: 001462  
Lot No.: 01834176  
Mobile Phase: 50:50 v/v Acetonitrile/0.10 M NH4OAc, pH 5.4  
Flow Rate: 1.00 mL/min  
Detection: UV, 254 nm  
Injection Volume: 5.0 µL  
Temperature: 30 °C  
Storage Solution: Acetonitrile

No.	Peak Name	Ret.Time (min)	Asymmetry (EP)	Efficiency (EP)	Concentration (µg/mL)
1	Cytosine	2.6	1.11	12106	100
2	Naphthalene	8.7	0.98	19685	100

#### QA Results:

Analyte	Parameter	Specification	Results
Naphthalene	Efficiency	>=16,200	Passed
Naphthalene	Asymmetry	0.95-1.32	Passed
Naphthalene	Retention Time	8.0-9.6	Passed
Naphthalene	Pressure	<=1980	1352

Production Reference:  
Database: Silica  
Directory: Silica/Silica\_1  
Sequence: 1533273\_068352\_AM  
Sample No: 12  
6.80 DU134 Bulk 3981 (22475) (Demo-Installation)  
Chromatop™ Thermo Fisher Scientific


## Acclaim Trinity P1

### Quality Assurance Report

Better Separations Through  
Better Chemistry

**Acclaim™ Trinity™ P1**  
3µm (3.0 x 150 mm)  
Product No. 075563

Date: 26-Jun-19 09:08  
Serial No.: 001444  
Lot No.: 01705135  
Mobile Phase: 50:50 v/v Acetonitrile/30 mM (total) ammonium acetate, pH 5  
Flow Rate: 0.50 mL/min  
Detection: UV, 230 nm  
Injection Volume: 2.0 µL  
Temperature: 30 °C  
Storage Solution: Mobile Phase

No.	Peak Name	Ret.Time (min)	Asymmetry (EP)	Efficiency (EP)	Concentration
1	1,1-Diethyl bisulfide	3.4	1.17	11972	0.02 mg/mL
2	Triphenylene	7.1	1.36	6154	0.02 mg/mL
3	Nitraz	11.2	1.08	12907	0.10 mg/mL

(1) Efficiency (Plates/Column)

#### QA Results:

Analyte	Parameter	Specification	Results
Nitraz	Efficiency	>=10,000	Passed
Nitraz	Asymmetry	0.95-1.28	Passed
Nitraz	Retention Time	10.1-12.9	Passed
1,1-Diethyl bisulfide	Asymmetry	<=1.7	Passed
1,1-Diethyl bisulfide	Pressure	<=2970	1434

Production Reference:  
Database: Silica  
Directory: Silica/Silica\_1  
Sequence: 1521599\_075563\_KM  
Sample No: 1  
6.80 SRI 4 Bulk 4521 (23099) (Demo-Installation)  
065167-02 (QAR)  
Chromatop™ Thermo Fisher Scientific

thermoscientific.com/dionex  
© 2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.  
XX2155-EN 02155 031298-11

Thermo Fisher Scientific  
12291 River Way  
P.O. Box 3003  
Sunnyvale, CA 94088-3003  
(408) 737-0700

**thermo**  
scientific

thermoscientific.com/dionex  
© 2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.  
XX2155-EN 02155 031298-11

Thermo Fisher Scientific  
12291 River Way  
P.O. Box 3003  
Sunnyvale, CA 94088-3003  
(408) 737-0700

**thermo**  
scientific

thermoscientific.com/dionex  
© 2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.  
XX2155-EN 02155 031298-11

Thermo Fisher Scientific  
12291 River Way  
P.O. Box 3003  
Sunnyvale, CA 94088-3003  
(408) 737-0700

**Thermo**  
SCIENTIFIC  
A Thermo Fisher Scientific Brand

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### Quality Assurance Report (QAR)

	出厂溶剂	活化方法	短期储存	长期储存
Column	Storage Solution	Column Flushing	For short-term storage (< 3 days)	For long-term storage
Acclaim Mixed-Mode WAX-1	50:50 v/v Acetonitrile / 50 mM Phosphate buffer, pH 6.0	50:50 Acetonitrile/20 mM Ammonium acetate, pH 5.0 at least 40 column volumes	Mobile phase at room temperature	Acetonitrile at room temperature
Acclaim Mixed-Mode WCX-1	Acetonitrile			80:20 Acetonitrile/10 mM Ammonium acetate, pH 5.0 at room temperature
Acclaim Mixed-Mode HILIC-1	60:40 v/v Acetonitrile / 100 mM Ammonium acetate, pH 5.0			
Acclaim Trinity P1	50:50 v/v Acetonitrile / 30 mM (Total) Ammonium acetate, pH 5.0			
Acclaim Trinity P2	60:40 v/v Acetonitrile / 100 mM Ammonium formate, pH 3.65			
Acclaim Trinity Q1	75:25 v/v Acetonitrile/100 mM (Total) Ammonium acetate, pH 5.0			
Acclaim Surfactant Acclaim Surfactant Plus	Acetonitrile			

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ThermoFisher  
SCIENTIFIC

Development guide | Hydrophobic columns

ThermoFisher  
SCIENTIFIC

## LC columns

### Column care guide and general method development information for Thermo Scientific hydrophobic columns

Applies to columns bonded with phases such as C18, C8 C4, C30 and polar embedded C18s, polar endcapped C18s and other proprietary hydrophobic phases

#### Before you get started

Manuals, specification sheets or technical guides for your column might be available to download from [thermofisher.com](http://thermofisher.com). Type the P/N or product name in the search box. Helpful literature is near the bottom of the product page. Some columns include a Quick-Start Guide in the box and/or a yellow caution tag on the column. Please read these before using the column.

Always start by investigating the Certificate of Analysis (CoA) or Quality Assurance Report (QAR) accompanying your column. This document includes a lot of valuable information. For instance, investigate what solvent the column is shipped in. If the column is filled with something incompatible with your mobile phase, flush it out with a mutually compatible intermediate solvent. Some detectors such as charged aerosol and mass spectrometers are highly sensitive to column bleed. Condition the column before connecting it to the detector.

You should always strive to reproduce the chromatogram in your CoA or QAR when you receive the column into your lab. This way you can assure that the column is operating correctly when you start your method, and if you routinely repeat the column's CoA or QAR, you can notice column degradation early on and implement preventative measures if needed.

For UH-PLC columns operating at high pressure > 400 bar, it can take the column 20–30 minutes of extra time to come to thermal steady-state after the column oven is ready. Continue equilibration until the pressure and detector baselines are stable.

Always check for leaks before use.

#### Operational limits

Respect the limits for pressure, pH, temperature and solvent compatibility. The product manual, specification sheet or technical guide is the best reference for operational limits. If there is not a manual, see the online [catalog](http://catalog) or product web page on [thermofisher.com](http://thermofisher.com)

Operating near the extremes of the pH or temperature limits can reduce column life and increase column bleed.

thermo scientific

Development guide | Phenyl columns

ThermoFisher  
SCIENTIFIC

## LC columns

### Column care guide and general method development information for Thermo Scientific phenyl columns

Applies to columns bonded with aromatic phases such as Phenyl, Phenyl X, Phenyl Hexyl, PFP and Biphenyl

#### Before you get started

Manuals, specification sheets or technical guides for your column might be available to download from [thermofisher.com](http://thermofisher.com). Type the P/N or product name in the search box. Helpful literature is near the bottom of the product page. Some columns include a Quick-Start Guide in the box and/or a yellow caution tag on the column. Please read these before using the column.

Always start by investigating the Certificate of Analysis (CoA) or Quality Assurance Report (QAR) accompanying your column. This document includes a lot of valuable information. For instance, investigate what solvent the column is shipped in. If the column is filled with something incompatible with your mobile phase, flush it out with a mutually compatible intermediate solvent. Some detectors such as charged aerosol and mass spectrometers are highly sensitive to column bleed. Condition the column before connecting it to the detector.

You should always strive to reproduce the chromatogram in your CoA or QAR when you receive the column into your lab. This way you can assure that the column is operating correctly when you start your method, and if you routinely repeat the column's CoA or QAR, you can notice column degradation early on and implement preventative measures if needed.

For UH-PLC columns operating at high pressure > 400 bar, it can take the column 20–30 minutes of extra time to come to thermal steady-state after the column oven is ready. Continue equilibration until the pressure and detector baselines are stable.

Always check for leaks before use.

#### Operational limits

Respect the limits for pressure, pH, temperature and solvent compatibility. The product manual, specification sheet or technical guide is the best reference for operational limits. If there is not a manual, see the online [catalog](http://catalog) or product web page on [thermofisher.com](http://thermofisher.com)

Operating near the extremes of the pH or temperature limits can reduce column life and increase column bleed.

thermo scientific

Development guide | HILIC columns

ThermoFisher  
SCIENTIFIC

## LC columns

### Column care guide and general method development information for Thermo Scientific HILIC columns

Applies to columns bonded with phases such as HILIC, amino, cyano and bare silica

#### Before you get started

Manuals, specification sheets or technical guides for your column might be available to download from [thermofisher.com](http://thermofisher.com). Type the P/N or product name in the search box. Helpful literature is near the bottom of the product page. Some columns include a Quick-Start Guide in the box and/or a yellow caution tag on the column. Please read these before using the column.

Always start by investigating the Certificate of Analysis (CoA) or Quality Assurance Report (QAR) accompanying your column. This document includes a lot of valuable information. For instance, investigate what solvent the column is shipped in. If the column is filled with something incompatible with your mobile phase, flush it out with a mutually compatible intermediate solvent. Some detectors such as charged aerosol and mass spectrometers are highly sensitive to column bleed. Condition the column before connecting it to the detector.

You should always strive to reproduce the chromatogram in your CoA or QAR when you receive the column into your lab. This way you can assure that the column is operating correctly when you start your method, and if you routinely repeat the column's CoA or QAR, you can notice column degradation early on and implement preventative measures if needed.

For UH-PLC columns operating at high pressure > 400 bar, it can take the column 20–30 minutes of extra time to come to thermal steady-state after the column oven is ready. Continue equilibration until the pressure and detector baselines are stable.

Always check for leaks before use.

#### Operational limits

Respect the limits for pressure, pH, temperature and solvent compatibility. The product manual, specification sheet or technical guide is the best reference for operational limits. If there is not a manual, see the online [catalog](http://catalog) or product web page on [thermofisher.com](http://thermofisher.com)

Operating near the extremes of the pH or temperature limits can reduce column life and increase column bleed.

thermo scientific

Hydrophobic LC Columns Care Guide-CN  
Phenyl LC Columns Care Guide-CN  
HILIC LC Columns Care Guide-CN

Mixed Mode LC Columns Care Guide-CN  
Hypercarb LC Columns Care Guide-CN  
HyperREZ LC Columns Care Guide-CN

# Hypersil GOLD 柱效测试报告

ThermoFisher  
SCIENTIFIC

thermoscientific

Part Number: 25005-254630  
Column: Hypersil GOLD™  
Serial Number: 20031942  
Lot Number: 16961  
Column Dimensions: 250 mm x 4.6 mm


PN / SN / LN

## Chromatographic Parameters

Mobile Phase: 60/40 Acetonitrile/Water  
Flow Rate: 1.25 mL/min  
Sample Volume: 2.5 µL  
Wavelength: UV @ 254 nm  
Particle Size: 5 µm  
Pore Size: 175 Å  
Temperature: Ambient  
Column Storage: Mobile Phase  
Column Back Pressure: 1650 psi

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
$$\epsilon_T = 2.835 / 4.135 \times 100\% = 68.3\%$$

The void volume of the column is the volume that is not taken up by the stationary phase (approximately 68% of the column volume):

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





### Chromatographic Parameters

Mobile Phase: 50/50 Acetonitrile/Water  
Flow Rate: 1.5 mL/min  
Temperature (°C): 30  
Sample Volume: 1.0 µL  
Wavelength: UV @ 254 nm

### Chromatographic Test

	Specification		Result
	Minimum	Maximum	
Retention Time- min (o-Xylene)	6.9	9.6	8.45
Capacity Factor (o-Xylene)	4.5	6.3	5.55
Efficiency- N/m (o-Xylene)	150000	220000	181972
USP Tailing Factor (o-Xylene)	0.9	1.2	1.05
Column Back Pressure* - psi	2500	4000	2938


\*These values are for column testing in the factory. Actual upper pressure limit for the column is higher. Please check the Technical Guide for more details. To convert PSI to Bar divide by 14.5.

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